

simex M20Plus portable suction pump



Item No.	100730
Flow rate	22 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg)
Nominal Voltage (PSU)	In: AC 100-240 V ~ 50-60 Hz / 1.5 A Out: DC 12 V / 5 A
Degree of protection pursuant to IEC 60601-1	Type BF
Protection class pursuant to IEC 60601-1	II
Operating time	40 minutes ON; 10 minutes OFF
Sound emission	Without vacuum: approx. 66 dB(A)
Dimensions (basic devices H x W x D)	290 x 259 x 130 mm
Weight (basic device)	Approx. 3,1 kg
CE mark	CE0843