

simex M30D portable suction pump



Item No.	100406-3
Flow rate	26 l/min ± 2 l/min (low flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg)
Nominal Voltage (PSU)	In: AC 100 – 240 V ~ 50-60 Hz / 1.5 A Out: DC 12 V / 5 A
Degree of protection pursuant to IEC 60601-1	Type BF
Protection class pursuant to IEC 60601-1	II
Operating time	40 minutes ON, 10 minutes OFF
Sound emission	without vacuum: approx. 62 dB(A)
Dimensions (basic devices H x W x D)	290 x 259 x 130 mm
Weight (basic device)	Approx. 3,7 kg
CE mark	CE0843