

simex S20 portable suction pump



Item No.	100400-3
Flow rate (measuring point tube attachment)	18 l/min $\pm$ 2 l/min (low flow)
Vacuum	max. -90 kPa ($\sim$ -900 mbar; -675 mmHg)
mains voltage, mains frequency, power consumption, current consumption at 230 V	230 V AC~ 50 Hz/ 80 W/ 0.35 A
Degree of protection pursuant to IEC 60601-1	Type BF
Protection class pursuant to IEC 60601-1	II
Operating time	continuous
Sound emission	without vacuum: approx. 55 dB(A)
Dimensions (basic devices H x W x D)	290 x 259 x 130 mm
Weight (basic device)	Approx. 3,4 kg
CE mark	CE0843