

simex S20K portable suction pump



| | |
|--|--|
| Item No. | 100402-3 |
| Flow rate (measuring point tube attachment) | 18 l/min ± 2 l/min (low flow) |
| Vacuum | max. -90 kPa (~ -900 mbar; -675 mmHg) |
| mains voltage, mains frequency, power consumption, current consumption at 230 V | 230 V AC~ 50 Hz/ 156 W/ 0.72 A |
| Degree of protection pursuant to IEC 60601-1 | Type BF |
| Protection class pursuant to IEC 60601-1 | II |
| Operating time | continuous |
| Sound emission | without vacuum: approx. 66 dB(A) |
| Dimensions (basic devices H x W x D) | 290 x 259 x 130 mm |
| Weight (basic device) | Approx. 3,5 kg |
| CE mark | CE0843 |