

## *simex S30 portable suction pump*



<b>Item No.</b>	<b>100401-3</b>
<b>Flow rate (measuring point tube attachment)</b>	<b>28 l/min ± 2 l/min (high flow)</b>
<b>Vacuum</b>	<b>max. -90 kPa (~ -900 mbar; -675 mmHg)</b>
<b>mains voltage, mains frequency, power consumption, current consumption at 230 V</b>	<b>230 V AC~ 50 Hz/ 80 W/ 0.35 A</b>
<b>Degree of protection pursuant to IEC 60601-1</b>	<b>Type BF</b>
<b>Protection class pursuant to IEC 60601-1</b>	<b>II</b>
<b>Operating time</b>	<b>45 minutes ON, 10 minutes OFF</b>
<b>Sound emission</b>	<b>without vacuum: approx. 57 dB(A)</b>
<b>Dimensions (basic devices H x W x D)</b>	<b>290 x 259 x 130 mm</b>
<b>Weight (basic device)</b>	<b>Approx. 4,1 kg</b>
<b>CE mark</b>	<b>CE0843</b>