

# SIM<sup>EX</sup> 200 Negative Pressure Wound Therapy



<b>Item No.</b>	<b>100139-3</b>	
<b>Flow rate</b>	<b>max. 8 l/min</b>	
<b>(Negative) Pressure</b>	<b>Max. -200mmHg</b> <b>Conversion factor: 1kPA~7,5mmHg</b>	
<b>Power requirement</b>	<b>100 – 240V, 50/60Hz, 15W</b>	
<b>Current consumption</b>	<b>1,25A</b>	
<b>Measurements (H x W x D)</b>	<b>165 x 220 x 90 mm</b>	
<b>Weight (basic unit)</b>	<b>Approx. 1,2 kg</b>	
<b>Running time</b>	<b>Mains</b>	<b>continuous operation</b>
	<b>Battery</b>	<b>Approx. 24-48 hours (depending on use)</b>
<b>Sound emission</b>	<b>35 dB (A)</b>	
<b>Degree of protection acc. IEC 60601-1</b>	<b>Typ BF</b>	
<b>Protection class acc. IEC 60601-1</b>	<b>II</b>	
<b>CE mark</b>	<b>CE0843</b>	