

SIM^{EX} 300 Negative Pressure Wound Therapy



Item No.		100399-3
Flow rate		max. 8 l/min
(Negative) Pressure		Max. -200mmHg Conversion factor: 1kPA~7,5mmHg
Power requirement		100 – 240V, 50/60Hz, 15W
Current consumption		1,25A
Measurements (H x W x D)		290 x 259 (+100 canister) x 130 mm
Weight (basic unit)		Approx. 2,2 kg
Running time	Mains	continuous operation
	Battery	Approx. 24-48 hours (depending on use)
Sound emission		35 dB (A)
Degree of protection acc. IEC 60601-1		Typ BF
Protection class acc. IEC 60601-1		II
CE mark		CE0843