

## *simex M30Plus portable suction system*



<b>Item No.</b>	<b>100731</b>
<b>Flow rate</b>	<b>30 l/min ± 2 l/min (high flow)</b>
<b>Vacuum</b>	<b>max. -90 kPa (~ -900 mbar; -675 mmHg)</b>
<b>Nominal Voltage (PSU)</b>	<b>In: AC 100 – 240 V~ 50-60 Hz / 1.5 A Out: DC 12 V / 5 A</b>
<b>Protection class acc. IEC 60601-1</b>	<b>Type BF / II</b>
<b>Ingress Protection (IP) acc. IEC 60529</b>	<b>IP22</b>
<b>Operating time</b>	<b>40 minutes ON, 10 minutes OFF</b>
<b>Battery, rechargeable</b>	<b>10.8 V; 5800 mAh (lithium ion battery)</b>
<b>Sound emission</b>	<b>without vacuum: approx. 66 dB(A)</b>
<b>Dimensions (basic devices H x W x D)</b>	<b>290 x 259 x 130 mm</b>
<b>Weight (basic device)</b>	<b>Approx. 3,6 kg</b>
<b>CE mark</b>	<b>CE0123</b>