

simex S20K portable suction system



Item No.	100402-3
Flow rate (measuring point tube attachment)	18 l/min ± 2 l/min (low flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg)
mains voltage, mains frequency, power consumption, current consumption at 230 V	230 V AC~ 50 Hz/ 156 W/ 0.72 A
Protection class acc. IEC 60601-1	Type BF / II
Ingress Protection (IP) acc. 60529	IP22
Operating time	continuous
Sound emission	without vacuum: approx. 66 dB(A)
Dimensions (basic devices H x W x D)	290 x 259 x 130 mm
Weight (basic device)	Approx. 3,5 kg
CE mark	CE0123