

simex S30 portable suction system



Item No.	100401-3
Flow rate (measuring point tube attachment)	28 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg)
mains voltage, mains frequency, power consumption, current consumption at 230 V	230 V AC~ 50 Hz/ 80 W/ 0.35 A
Protection class acc. to IEC 60601-1	Type BF / II
Ingress Protection (IP) acc. 60529	IP22
Operating time	45 minutes ON, 10 minutes OFF
Sound emission	without vacuum: approx. 57 dB(A)
Dimensions (basic devices H x W x D)	290 x 259 x 130 mm
Weight (basic device)	Approx. 4,1 kg
CE mark	CE0123