

simex cuff M automated subglottic aspiration system



| | | |
|--|---|---|
| Item No. | 100678-3 | |
| Flow rate | max. 8 l/min | |
| Pressure | -60 mbar to -300 mbar (in steps of 10 mbar) | |
| Containers | Disposable secretion container (250 ml / 750 ml) | |
| Nominal voltage of the power supply unit (In / Out) | 100 – 240V AC / 12V DC | |
| Maximum load current | 2,0 A | |
| Mains frequency (mains-powered) | 50 / 60 Hz | |
| Power consumption at 12 V | 24 W | |
| Rechargeable battery | 7,4 V; 4400 mAh (lithium-ion battery) | |
| Dimensions (H x W x D) | 165 x 220 x 90 mm | |
| Weight (basic device) | 1,2 kg | |
| Operating mode | Intermittent Aspiration | |
| Running time | Mains | continuous operation |
| | Battery | approx. 18 h, depending on the strain of the motor |
| Sound emission (Operation) | 35 dB (A) | |
| Protection class acc. IEC 60601-1 | Type BF / II | |
| Ingress Protection (IP) acc. IEC 60529 | IP33 | |
| CE mark | CE0123 | |