

simex cuff S automated subglottic aspiration system



Item No.	100679-3	
Flow rate	max. 8 l/min	
Pressure	-60 mbar to -300 mbar (in steps of 10 mbar)	
Container	Disposable secretion container system, 1000 ml	
Nominal voltage of the power supply unit (In / Out)	100 – 240V AC / 12V DC	
Maximum load current	2,0 A	
Mains frequency (mains-powered)	50 / 60 Hz	
Power consumption at 12 V	24 W	
Rechargeable battery	7,4 V; 4400 mAh (lithium-ion battery)	
Dimensions (H x W x D)	290 x 259 + 100 (container) x 130 mm	
Weight (basic device)	2,2 kg	
Operating mode	Intermittent Aspiration	
Running time	Mains	continuous operation
	Battery	approx. 18 h, depending on the strain of the motor
Sound emission (Operation)	35 dB (A)	
Protection class acc. IEC 60601-1	Type BF / II	
Ingress Protection (IP) acc. IEC 60529	IP33	
CE mark	CE0123	