

simex M20 portable suction system



Item No.	100403-3
Flowrate	18 l/min ± 2 l/min
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg)
Nominal Voltage (PSU)	In: AC 100 – 240 V ~ 50-60 Hz / 1.5 A Out: DC 12 V / 5 A
Protection class acc. IEC 60601-1	Type BF / II
Ingress Protection class acc. IEC 60529	IP22
Operating time	45 minutes ON, 10 minutes OFF
Battery, rechargeable	10.8 V; 5800 mAh (lithium ion battery)
Sound emission	Without vacuum: approx. 62 dB(A)
Dimension (basic device: H x W x D)	290 x 259 x 130 mm
Weight (basic device)	approx. 3.1 kg
CE mark	CE0123